Permanent and Provisional Abutments and prosthetic accessories

Conventional
Impression coping

Multi-Unit
Abutments

NeoBase
Abutments

Overdenture Equator
Abutments

Prepable Titanium
Abutments

English (English)	Instructions for use Permanent and Provisional Abutments and prosthetic accessories.....	2
Српски (Serbian)	Uputstvo za upotrebu Trajne i privremene suprastrukture i protetski pribor.....	9
Dansk (Danish)	Brugsanvisning Permanente og midlertidige brotænder og protesetilbehør.....	16
Deutsch (German)	Gebrauchsanweisung Definitive und provisorische Abutments und prothetisches Zubehör.....	23
Ελληνικά (Greek)	Οδηγίες χρήσης Μόνιμα και προσωρινά κολοβώματα και προσθετικά εξαρτήματα.....	30
Español (Spanish)	Instrucciones de uso Pilares provisionales y definitivos y accesorios protésicos.....	37
Français (French)	Notice d'utilisation Piliers permanents et provisoires et accessoires prothétiques.....	44
Italiano (Italian)	Istruzioni per l'uso Pilastrì permanenti e provvisori e accessori protesici.....	51
Magyar (Hungarian)	Használati útmutató Állandó és ideiglenes felépítmények és protézistartozékok.....	58
Nederlands (Dutch)	Gebruiksaanwijzing Permanente en voorlopige abutments en prothetische accessoires.....	65
Norsk (Norwegian)	Bruksanvisning Permanente og provisoriske distanser og protesetilbehør.....	72
Polski (Polish)	Instrukcja obsługi łączniki stałe i tymczasowe oraz akcesoria protetyczne.....	79
Português (Portuguese)	Instruções de utilização Pilares provisórios e definitivos e acessórios protéticos.....	86
Svenska (Swedish)	Bruksanvisning Permanenta och temporära distanser och protetiska tillbehör.....	93

Permanent and Provisional Abutments and prosthetic accessories

Conventional
Impression coping

Multi-Unit
Abutments

NeoBase
Abutments

Overdenture Equator
Abutments

Prepable Titanium
Abutments

Products (Product Categories)

Impressions: Conventional Impression coping · Equator Impression coping · ScanPeg · TiBase ScanPost · Scan Body

Temporary abutments: Esthetic Tissue Former (Provisional Abutment) · Provisional Titanium Abutment

Final abutments: Abutment Blanks · Access Abutments · CoCr Abutments · Gold Abutments · Multi-Unit Abutments · NeoBase Abutments · Neoss TiBase Abutments · Overdenture Ball Abutments · Overdenture Equator Abutments · Prepable Titanium Abutments · Prepable Zirconia Abutments · Titanium NeoLink Abutments

Abutment/prosthetic screw: Abutment/Prosthetic screw iGO · Abutment/Prosthetic Screw Neo

Burnable abutment · Burnable coping: Burnout Abutment · Burnout Coping

Laboratory: Laboratory screw iGO · Laboratory screw Neo · Replica

Description

The Neoss implant system is a complete system of endosseous dental implants with corresponding instruments, prosthetic products and accessories.

For additional information on the specific Neoss Implant System product you are using please consult the Neoss Implant System Clinical Handbook (10501) and translated versions.

Permanent and provisional abutments and prosthetic accessories are designed for prosthetic restorations and impression taking on Neoss dental implants and abutments. These abutments and prosthetic accessories are available in a variety of shapes, sizes and materials to fit individual patient needs and different production techniques.

For a specific product description, article number and dimensions please consult the product label and the Neoss Product Catalogue.

Intended use

Neoss abutments are provisionally or permanently connected to Neoss implants or Multi-Unit or Access Abutments and used as the base onto which a fixed (non-removable) prosthetic restoration is fitted for patient rehabilitation. They also support and shape the soft tissue healing.

Abutment screws are used to fixate the abutment to the implant. Prosthetic screws are used to fixate a fixed prosthetic restoration to the Multi-Unit abutment or Access abutment.

Neoss Overdenture Abutments are intended as support for a removable permanent denture, existing or new. Overdenture attachments are embedded into the denture to keep denture in place.

The purpose of an impression is to accurately transfer the position of the implant or the abutment, in relation to the patient specific anatomical situation to a model for accurate manufacturing of the prosthesis. The ScanPeg is a scan body for intraoral scanning momentarily fitted in the screw access hole of the Esthetic Healing Abutment to enable digital acquisition of the implant position in relation to the adjacent teeth and soft tissue.

Non-clinical Laboratory and Burnout products are used in the dental laboratory to manufacture prosthetics that fit correctly onto Neoss implants or Multi-Unit or Access abutments.

Abutments and Impression products (except for TiBase Impression) are intended for single use for a single patient.

Laboratory screws and replicas are intended for single use.

Indication for use

Abutments and prosthetic accessories are used in conjunction with Neoss dental implants. Neoss implant system is intended to be used for treatment of edentulism by the replacement of one or more missing or failing teeth.

- Implants can be placed in fresh extraction sockets (immediate placement) as well as in healed sites.
- Implant can be placed using one-stage or two-stage surgical procedures.
- Immediate, early or delayed loading procedure can be applied.
- Fixed as well as removable prosthetics may be used.
- Prostheses may be single implant or splinted restorations.

Recognizing sufficient bone volume and stability and appropriate occlusal loads for the chosen procedure.

Patient target group

Adult population. Patients that have reached skeletal maturity and their jawbone is fully developed.

For all abutments, patients treated with Neoss dental implants.

Intended user and clinical setting

A licensed dental practitioner or physician will perform the clinical treatment using standardized, well-established surgical or prosthetic procedures in a hospital/clinic setting.

The fabrication of the restoration is performed in a dental laboratory setting by a licensed dental technician.

Seating of abutments and restorations in the patient's mouth take place in a dental office or in an environment of a comparable clinical standard by a licensed dental practitioner or physician.

Expected clinical benefits

The product is expected to fulfil its intended use in the specified indications, in its patient target group, during its expected lifetime.

Expected lifetime of the product in function

Permanent abutments, permanent abutment/prosthetic screws

The patient can, if not belonging to any groups with contraindications and if product is used as intended and for correct indications, expect the dental implant treatment to restore and maintain aesthetics and/or function for at least 25 years.

Provisional Abutments and Provisional Abutment/Prosthetic Screws

The patient can, if not belonging to any groups with contraindications and if the product is used as intended and for correct indications, expect the products to stay in place until replaced with abutments and restorations.

Overdenture Attachments

With repeated snapping in and out of the denture prosthesis, the Overdenture attachments eventually wear out over time, resulting in diminished snap retention and perceived looseness of the overdenture. When this happens, they are to be replaced with new attachments.

Non-implantable Single use products

Expected lifetime is one treatment in one patient.

For non-clinical single use products the expected lifetime is to work as intended to prepare one prosthetic restoration.

Reusable products

Can be reused as long as the integrity and performance of the products are maintained. The products must be inspected before each use for visible signs of damage, deformation, wear or corrosion. Products showing any signs of visible damage or loss of functional compatibility (for example fit of instrument in mating part) shall be discarded. Generally refer to instruction 14077 for assessment.

Contraindications

- Patients who are medically unfit for the medical procedure.
- Placement of implants in insufficient numbers or sizes to support biomechanical loads or undesirable positioning is contraindicated. Influencing factors might be narrow diameter implants, angulated abutments, posterior region, insufficient bone quality and quantity. Such an example is narrow diameter implants in combination with angulated abutments such as angulated Access Abutments in the posterior region.
- Treatment is contraindicated where the patient has a pre-existing allergy to products used for the treatment.

The Access abutments and multi abutments are contraindicated for single implant prosthetic cases.

Angulation correction of the prosthesis compared to the implant axis of more than 30° using a NeoBase or 20° using a TiBase, should be avoided.

Single implant restorations whose length exceeds a ratio of 1:1.25 in comparison to the length of the implant, should be avoided for NeoBase and TiBase.

Neoss Overdenture Abutments and attachments are contraindicated for single implant prosthetic cases and on implants with a greater divergence than 28° (50° with Housing for extended divergence) in conjunction with Equator abutments and 10° with Ball abutments.

Material

Permanent abutments are made of commercially pure titanium grade 4, titanium alloy (Ti6Al4V), Gold, Cobalt-chromium alloy or Zirconia. Equator Abutments have a titanium nitride (TiN) coating.

Healing abutments are made of PEEK (Polyetheretherketone).

Abutment screws are made of titanium alloy (Ti6Al4V) with a pure gold deposition.

Provisional abutments are made of titanium grade 4 or PMMA (Esthetic Tissue Formers).

Provisional screws are made of titanium alloy (Ti6Al4V).

Impression copings and posts are made of stainless steel or commercially pure titanium grade 4 if blue colour coded for NP components.

ScanPegs and Burnout products are made of PMMA.

Replicas are made in stainless steel.

Scan Body is made of titanium alloy (Ti6Al4V).

Sterility at delivery

Most abutments, screws, prosthetic accessories and non-clinical products are typically not provided sterile.

Products that are supplied sterile, such as Access abutments, are marked sterile on the packaging with a given expiry date. Sterile abutments are irradiated using Gamma.

Damaged packaging

If the packaging of a non-sterile product is damaged make sure to inspect the product visually and ensure that the product itself isn't damaged.

If the packaging of a sterile product is damaged upon delivery the product must be disposed.

Preparations before use

Preparation before use is required for Provisional abutments, Prepable abutments, Gold abutments, CoCr abutments and Abutment blanks according to Neoss Implant System Clinical Handbook (10501).

Sterile products, products for impression taking, replicas, screws or overdenture abutments and attachments shall not be adjusted before clinical use.

Cleaning, disinfection and sterilisation

All prosthetic and reusable products for invasive use including accessories delivered non-sterile must be unpacked, cleaned and if required sterilized before use as per instructions in 14077. This also applies for products that have been delivered sterile but adjusted or processed before clinical application.

The fabricator (dental technician) and the restoration must inform the dentist of the need to sterilize the abutment before clinical application (for US only).

Storage

Unused products must be stored in their original packaging. Sterilized bags and unused products must be stored in dry environment at room temperature and out of direct sunlight.

Procedure

For detailed information on procedures for a specific Neoss ProActive Implant System product, please consult Neoss Implant System Clinical Handbook (10501) which also includes information on procedures related to the dental laboratory.

Impression

The clinician either takes a conventional impression using an impression material with the impression coping attached to the implant or abutment, or a digital impression using a compatible digital impression coping such as ScanPeg or TiBase coping as per the Neoss Implant System Clinical Handbook (10501).

The attachment of the impression coping is made by hand tightening or 'push' in depending on the design.

Models

Based on the impression, the dental technician creates a model, physical or digital, with incorporated replicas representing the implants or abutments used for the impression.

Preparation before clinical placement

The following abutments or copings are adjusted for the clinical situation or incorporated into a restoration according to Neoss Implant System Clinical Handbook (10501) prior to or in conjunction with the clinical placement using a replica and provided laboratory screw to facilitate the attachment to a model during adjustment:

CoCr Abutments · Gold Abutments · NeoBase Abutments · Neoss TiBase Abutments · Prepable Titanium Abutments · Prepable Zirconia Abutments · Titanium NeoLink Abutments · Provisional Abutment Tissue Former · Provisional Titanium Abutment · Burnout Abutment · Burnout Coping · Abutment Blanks

Depending on workflow, different design and production techniques are used including digital design tools and milling options referenced in Neoss Implant System Clinical Handbook (10501).

Products attached to the replica are hand-tightened using the lab screw. Design and preparation work is conducted by a dental technician, or for certain solutions chair side by a clinician, following product specific design limitations such as post height, post angulation, marginal height, wall thickness, material processing and cementation. Provisional restorations used during the healing phase, should be designed out of occlusion. Make sure to protect the implant connection on the abutment during blasting, polishing or other procedures that can damage the connection. Surfaces exposed to the oral cavity or soft tissue are recommended to be smooth.

Copings or crowns, temporary or permanent, are made following standard procedures according to the material manufacturer's instructions.

The restoration is forwarded to the clinician together with the abutment screw or provisional screw, packaged to avoid any damages during transportation.

Placement of abutment or restoration

Non-sterile delivered dental restorations including abutments and copings should be cleaned, disinfected and when required sterilized as per the "Cleaning and disinfection" and "Sterilization" sections above prior to being clinically used.


Confirm that the appropriate Neoss abutment screw is used (applicable for all abutments except Overdenture abutments where the screw is incorporated into the abutment). The Neo screwdriver is used for all abutments except for abutments with angulated screw channels (NeoBase ASC and Abutment blanks ASC) where instead the iGO screwdriver is required.

After the removal of the cover screw, healing abutment or provisional restoration from the implant, clean and dry the implant connection.


Position the abutment or single tooth restoration in the implant while making sure that the retentive elements of the implant/abutment connection are properly aligned before the screw is tightened.

Tighten the abutment screw using the screwdriver in combination with the Neoss Ratchet set at the required torque listed below and on the packaging of the abutment or the screw. Note that Access angulated abutments requires a Neo screwdriver of 22 mm or longer.

Screw and Screwdriver compatibility

iGO Screws and iGO Screwdriver

Neoss Implant System Torque Recommendation (Ncm)

Provisional		Impression	Final restoration							
Esthetic Tissue Former	Provisional Ti Abutment	Impression Coping Screw	Neo Abutment Screw	iGO Abutment Screw		Multi-Unit and Access Abutment Straight & Angulated	Multi-Unit and Access Neo Prosthetic Screw	Ball Abutment including Access level	Equator Abutment	
				SP	NP				Implant level	Access level
20	20	10 Max	32*	32	20	32	20	20	32	20

* 35 Ncm optional for high-load cases

Warnings

The Neoss implant system should only be used by a licensed dental practitioner or physician who have had the appropriate education and training.

Implement proper radiographic examination and planning to avoid vital anatomical structures (i.e. nerves, blood vessels, teeth or other sensitive structures) during implant site preparation and implant insertion.

Because of the small size of the products, care must be taken that they are not swallowed or aspirated by the patient.

The implant to the abutment connection is essential to the mechanical stability of the dental implant system. The abutment connection is not to be modified.

Any change to this abutment connection will characterize your facility as a medical device manufacturer subject to FDA registration, fees, regulation and restrictions (for US only).

Precautions

- Neoss implant system products should only be used together with original products and instruments.
- Routine implant treatment is not recommended in pediatric patients with uncompleted maxillary and mandibular growth.
- Neoss implants, cover screws, abutments and abutment screws must be used solely on one patient. Single use products should not be reused due to risks of product contamination, patient/user infection and/or failure of the product to perform as intended.
- All multiple use products must be maintained in good condition and inspected before use to avoid harm to patient and damage to products. See Document 14077 (Cleaning and maintenance of Neoss Reusable Products) for inspection guidelines.
- Careful clinical and radiographic examination to assess the overall medical status of the patient is performed before treatment. Examples of risk factors for dental implant treatment are (but not limited to):
 - Insufficient bone quantity or quality
 - Local or systemic infections or inflammation
 - Compromised general and local health
 - Poor oral hygiene
 - Smoking – Alcohol or drug abuse
 - Disorders, medications or therapies affecting bone and wound healing (i.e. diabetes, chemotherapy, bisphosphonates)
 - History of or ongoing therapeutic radiation in the area
 - Uncontrolled bleeding disorders

- Uncontrolled parafunctional habits or unfavorable jaw relationships
- Poor patient compliance
- Assess vertical and horizontal intraoral space prior to treatment to assure unobstructed use of all components and instruments.
- Biomechanical principles should be recognized when choosing loading protocol as well as for the design and construction of the abutment and prosthesis to minimize unfavorable forces. Failure to do so can lead to excessive bending force and fatigue failure of the implant or abutment components.
- When possible, modifications to prosthetic reconstructions should take place extraorally.
- Waste materials should be handled according to the established procedures for hazardous material at the hospital/clinic/laboratory.
- In the event of malfunction of the product or changes in its performance of the product, the patient should contact the dentist for assessment.

Adverse effects

Adverse effects related to dental implant treatment includes but are not limited to:

Intervention related trauma · Damage to adjacent teeth · Fracture of the jaw bone · Temporary or permanent nerve injury (dysesthesia or paresthesia) · Unintended sinus membrane perforation · Post-operative bleeding · Hematoma · Aspirated or swallowed products · Allergic reaction · Malpositioned implant · Postoperative discomfort/swelling/pain · Bone necrosis · Infections · Inflammatory reactions · Pain · Implant loss/failure/fracture · Prosthetic product failure/fracture · Excessive bone loss · Excessive gingival pressure/tension · Peri-implantitis · Phonetic difficulties · Esthetic problems

Magnetic resonance imaging (MRI)

MRI Safety Information	 MR Conditional
A patient with a system out of bridges, bars, and implants from Neoss AB may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Static Magnetic Field Strength (B0) [T]	1.5T or 3.0T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Maximum Whole-Body SAR	1 W/kg 1.5T and 3T without any position restriction 2 W/kg 1.5T and 3T nearest part of implant \geq 25 cm out of isocenter feet-head direction
Maximum Head SAR	1.6 W/kg
Scan Duration and Wait Time	Continuous scanning for 60 minutes
MR Image Artifact	The presence of this device may produce an image artifact
<i>Additional information about Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Tests following guidelines in ASTM F2119-07 (2013) showed artifact widths of \leq5 mm from the surface of the device for the gradient echo sequence and \leq4 mm from the surface of the device for the spin echo sequence.</i>	

For the USA only

Caution: Federal (USA) law restricts the sale of this product to or on the order of a licensed physician or dentist.

Additional information


Serious incidents

If any serious incident occurs in relation to the product, the user and/or patient should report to the manufacturer (<https://www.neoss.com>) and the competent authority of the state in which the user and/or patient is established.

Safe disposal of the product

Handling of hazardous material according to established procedures at the hospital/clinic. The disposal of the product shall be performed in an environmentally sustainable manner according to local regulations. If the product is contaminated with human blood, tissue residues or other human secretions the product shall be disposed in appropriate containers for this.

Symbols used on packaging

NON-STERILE 	CAUTION 	KEEP DRY 	USE BY/ EXPIRY DATE 	NARROW PLATFORM 	STANDARD PLATFORM 
UNIQUE DEVICE IDENTIFIER 	LOT/BATCH NUMBER 	CONSULT INSTRUCTIONS FOR USE 	DATE OF MANUFACTURE 	MANUFACTURER 	
MEDICAL DEVICE 	CATALOGUE NUMBER 	STERILIZED USING IRRADIATION 	STERILIZED USING ETHYLENE OXIDE 	KEEP AWAY FROM SUNLIGHT 	
DO NOT USE IF PACKAGE IS DAMAGED 	DO NOT RE-USE (single use only) 	CAUTION: Federal (USA) law restricts the sale of this device to or on the order of a licensed physician or dentist Rx only	CE MARK AND NOTIFIED BODY NUMBER CE 2797		
SINGLE STERILE BARRIER SYSTEM 	SINGLE STERILE BARRIER SYSTEM with protective packaging outside 	SINGLE STERILE BARRIER SYSTEM with protective packaging inside 	CONTAINS HAZARDOUS SUBSTANCES (Co) 		

Trajne i privremene suprastrukture i protetski pribor

Konvencionalna
osa otiska

Multi-Unit
suprastrukture

NeoBase
suprastrukture

Equator suprastrukture
za proteze

Prepable titanijumske
suprastrukture

Proizvodi (kategorije proizvoda)

Otisci: Konvencionalna osa otiska · Osa otiska Equator · ScanPeg · TiBase ScanPost · Skeniranje tela

Privremene suprastrukture: Esthetic Tissue Former (Provisional Abutment) · Provisional Titanium Abutment

Trajne suprastrukture: Upornice · Access suprastrukture · CoCr suprastrukture · Gold suprastrukture · Multi-Unit suprastrukture · NeoBase suprastrukture · Neoss TiBase suprastrukture · Ball suprastrukture za proteze · Equator suprastrukture za proteze · Prepable titanijumske suprastrukture · Prepable cirkonske suprastrukture · Titanijumske NeoLink suprastrukture

Vijak suprastrukture/protetski vijak: Vijak suprastrukture/protetski vijak iGO · Vijak suprastrukture/protetski vijak Neo

„Burnable“ suprastruktura · „Burnable“ kapica: „Burnout“ suprastruktura · „Burnout“ kapica

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Permanenta distanser: Abutment Blanks · Access Abutments · CoCr Abutments · Gold Abutments · Multi-Unit Abutments · NeoBase Abutments · Neoss TiBase Abutments · Overdenture Ball Abutments · Overdenture Equator Abutments · Prepable Titanium Abutments · Prepable Zirconia Abutments · Titanium NeoLink Abutments

Distans/proteesskruv: Abutment/Prosthetic screw iGO · Abutment/Prosthetic Screw Neo

Burnable distans · Burnable coping: Burnout Abutment · Burnout Coping

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