Drill Sleeve Ø2.2

This Instructions For Use (IFU) is for Neoss Drill Sleeve Ø2.2

Description

A drilling guide can be produced with an incorporated drill sleeve. Different production techniques can be used such as vacuum forming and conventional stent forming techniques or it can be produced with milling or printing on the basis of 3D data.

Indications

Neoss Drill Sleeves are used for insertion into planning and surgical drilling guides. Drill sleeves are indicated for pre-implant planning and initial drilling to mark the center point of the implant position.

Procedure – Guide manufacturing

Based on input from the clinician a drill guide is either made on a stone model or designed in a software and milled or printed.

Manual procedure using a stone model:

The guide structure is created by vacuum forming or other conventional stent forming techniques in a dental laboratory environment. Holes with same or similar outer dimensions as the drill sleeve are drilled in desired position and angulation. Compatible profile drills are recommended. The drill sleeves are then positioned either by gluing or press fit. The position of the sleeve is defined on the basis of anatomical and prosthetic requirements.

Based on digital processing and a scanned model:

The geometries of the drill sleeve are available in libraries of selected 3D implant planning and guide design software programs. Please contact your local Neoss representative for available options.

Once the guide is designed in the planning and guide software the drill sleeves can be chosen from different libraries. The position of the sleeve is defined on the basis of anatomical and prosthetic needs. The CNC milling machine or the 3D printer creates the holes for the sleeves in the guide. The sleeve is positioned into the guide according to the guide manufacturer's recommendation and fastened accordingly through press fit or gluing.

Procedure – Clinical application


Refer to cleaning and clinical placement guidance related to the specific guide solution used. With the guide in place the Ø2.2 mm Neoss drill is chosen and inserted through the guide. Once in place the drilling can start. As the guide is only for marking the entry point and initial directional guidance, the drill should be inserted to a shorter depth than the length of the implant.

Please note that the drill sleeve is not intended to control the depth, hence there is no drill stop function.

For additional information on Neoss drills and detailed information on the other specific Neoss Implant System products you are using, please consult general Instructions For Use (10538) and Neoss Implant System Guidelines (10501).

Sterilization

Neoss Drill Sleeves are provided non-sterile and can be sterilized by autoclaving. Please refer to IFU Cleaning and Maintenance of Neoss Reusable Products (14077).

Neoss products can withstand temperatures up to 150°C unless otherwise stated. This process only applies to the drill sleeves, please refer to material manufacturer for cleaning, disinfection and sterilization instructions for the guide material used.

The instruments and accessories are reusable devices that can be reused as long as the integrity and performance of the device are maintained. The devices must be inspected before each use for visible signs of damage, deformation, wear or corrosion. Devices showing any signs of visible damage or loss of functional compatibility (for example fit of instrument with mating components, lifting function etc) shall be discarded.

Storage

Sterilized bags and unused components must be stored in dry environment, at room temperature and out of direct sunlight. Unused components must be stored in their original packaging.

Contraindications

Treatment is contraindicated where the patient has a preexisting allergy to the used parts. Use of drill sleeves with other than Neoss drills is contraindicated. There are no other specific contraindications other than what is related to the general implant treatment.

Disposal

Handling of hazardous material according to established procedures at the hospital/clinic.

Cautions

- Federal (USA) law restricts the sale of this device to or on the order of a licensed physician or dentist.
- Drill Sleeves must be used solely on one patient. Single use devices should not be reused due to risks of product contamination, patient/user infection and/or failure of the device to perform as intended.
- Care must be taken that the Drill Sleeves are not swallowed or aspirated by the patient.
- In the event of malfunction of the device or changes in its performance of the device, the patient should contact the dentist for assessment.

General packaging symbols

USE BY/EXPIRY DATE 	CATALOGUE NUMBER 	LOT/BATCH NUMBER 	DO NOT RE-USE (Single use only) 	TEMPERATURE LIMIT 
KEEP AWAY FROM SUNLIGHT 	MANUFACTURER 	DATE OF MANUFACTURE 	STERILIZED USING ETHYLENE OXIDE 	STERILE BY IRRADIATION (Contents of inner package sterile) 
NON-STERILE 	DO NOT USE IF PACKAGE IS DAMAGED 	MEDICAL DEVICE 	CONSULT INSTRUCTIONS FOR USE (Also available on www.neoss.com/IFU) 	CAUTION: Federal (USA) law restricts the sale of this device to or on the order of a licensed physician or dentist Rx only

 **Neoss AB**
Arvid Wallgrens backe 20
413 46 Göteborg
Sweden
T +46 (0)31 88 12 80
W www.neoss.com

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