EN

Applications


The torque ratchet is designed for controlled manual insertion of implants and tightening abutment screws under a defined torque. The appropriate instrument (i.e. Manual Handle, Wrench Adapter or Manual Screwdriver) is inserted and carried by the ratchet head. The removal of the instrument becomes easier if the pin (5) is drawn away from the instrument. The preset torque is set using the adjusting nut (4).

Maintenance of the Torque Ratchet

In addition to the cleaning, disinfection, sterilization and storage instructions described in Cleaning and maintenance of Neoss Reusable Products (14077) please note the following:

- After use and in preparation for pre-cleaning, dismantle the ratchet – this does not require any tools.
- After cleaning and rinsing, all components are dried immediately.
- Prior sterilization, the marked areas with  should be slightly moistened with handpiece maintenance oil.
- Reassemble the ratchet parts in a relaxed position (setting about 10 Ncm). The labelling IN on the ratchet head (1) and scale (6) face the same direction.
- The lifetime of the ratchet is primarily dependent on care and not the number of sterilization cycles. The device can be reused as long as the integrity and performance of the device is maintained.
- The ratchet shall be disposed of if, for example, the parts are not moving smoothly, are difficult to dismantle or show signs of discoloration.

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