

Patient Information Leaflet NeoGen PTFE membrane

English (English)	Patient Information Leaflet NeoGen PTFE membrane	2
عربي (Arabic)	نشرة معلومات المريض للتعرف على غشاء البولي رباعي فلورو إيثيلين (PTFE) من NeoGen	4
Dansk (Danish)	Brochure med patientoplysninger til NeoGen PTFE-membran	6
Deutsch (German)	Beipackzettel NeoGen PTFE-Membran	8
Ελληνικά (Greek)	Ενημερωτικό φυλλάδιο ασθενούς για τη μεμβράνη PTFE της NeoGen	10
Español (Spanish)	Folleto de información para el paciente de la membrana de PTFE NeoGen	12
Français (French)	Notice d'information destinée aux patients sur la membrane en PTFE NeoGen	14
Italiano (Italian)	Foglietto illustrativo per il paziente – Membrana NeoGen in PTFE	16
Magyar (Hungarian)	NeoGen PTFE membrán betegtájékoztatója	18
Nederlands (Dutch)	Patiëntenbijsluiter NeoGen PTFE-membraan	20
Norsk (Norwegian)	Pasientinformasjonsbrosjyre, NeoGen PTFE-membran	22
Polski (Polish)	Ulotka z informacjami dla pacjenta dot. membrany NeoGen z PTFE	24
Português (Portuguese)	Folheto Informativo da Membrana NeoGen PTFE	26
Svenska (Swedish)	Patientinformationsblad för NeoGen PTFE-membran	28

Patient Information Leaflet NeoGen PTFE membrane

Device description

The membrane is made of polytetrafluoroethylene (pure PTFE), which is a plastic material.

The membrane is strengthened with a titanium layer (medical grade commercially pure titanium, grade 2 (10 – 30 % weight)).

Intended purpose

The membranes should be used when there is a need to increase the amount of bone in the jaw to allow procedures, such as placing a dental implant.

The membranes are temporary and are removed when sufficient amount of bone has been created. The membranes are placed between the gum and the jawbone in a surgical procedure. The membranes create space for new bone to grow.

The membrane can be placed either before or at the same time as a dental implant.

The expected lifetime for the membrane is 6 months. The membrane is removed during a surgical procedure when sufficient bone has been created and treatment is complete.

When should the membrane be used?

Any patient where the amount of bone in the jaw needs to be increased (for example to place dental implants).

When should the membrane not be used?

- Patients who are medically unfit for the treatment.
- Patients with known allergy to material used for the treatment (PTFE and titanium).
- The product should not be placed in an area of active infection.

Risks, warnings, and precautions

Contact your healthcare professional if you believe that you are experiencing side effects related to the device, its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

Dental membranes have been safely used in dental surgeries for many years. Possible complications that can occur include:

- Mild discomfort/swelling/pain after the surgical procedure
- Inflammation or infection
- Membrane is exposed in the mouth
- Allergic reactions
- Numbness
- Bleeding
- Recession of the gum
- Abscess
- Loss of bone

- Difficulties to remove the membrane after completed treatment. Membrane material (fibers) is being left in the wound
- Implant loss
- Nausea and dizziness

Examples of risk factors for membrane treatment are (but not limited to):

- Infections or inflammation
- Compromised general and local health
- Poor oral hygiene
- Smoking
- Alcohol or drug abuse
- Disorders, medications, or therapies affecting bone and wound healing
- History of or ongoing radiation therapy in the area
- Uncontrolled bleeding disorders
- Poor patient compliance

MRI (Magnetic Resonance Imaging)

If you need to undergo an MRI scan, make sure to inform the healthcare provider that you have a NeoGen membrane, and make sure to present the Implant Card. A patient with a membrane may be safely scanned under the conditions stated in the Instruction for use.

Serious incidents

If any serious incident occurs in relation to the product, the user and/or patient should report to the manufacturer (<https://www.neoss.com>) and the competent authority of the state in which the user and/or patient is established.

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Neoss AB