

Software reference list

COMPANY	PLANNING						GUIDED				
	Software	Neoss ProActive Implants					Guide	"Fully" guided	"Partially" guided	Ø3.5 – 4.5 Implants**	Ø5.0 – 6.0 Implants** (+Neoss sleeve)
		Straight	Tapered	Edge	Ø6.0	Ø6.5					
360imaging	360dps	☑	☑	–	☑	–	360ips	☑	☑	☑	☑
3DDX	Osseoview	☑	☑	In progress	☑	In progress	3DDX Surgical Guide	☑	☑	☑	☑
3DIEMME	3Diagnosys	☑	☑	☑	☑	☑	RealGUIDE	☑	☑	☑	☑
3Shape	Implant Studio	☑	☑	☑	☑	☑	Implant Studio	☑	For review	☑	☑
Acteon	3Dapp	☑	☑	☑	☑	☑	N/A				
Anatomage	Invivio 5	☑	☑	In progress			Anatomage Guide	☑	☑	☑	In progress
Blue Sky Bio	Blue Sky Plan	☑	☑	☑	☑	☑		–	–	–	–
Carestream	CS 3D Imaging	☑	☑	☑	☑	☑	N/A				
Centaur	Mediasuite 3D	☑	☑	–	☑	–	N/A				
Dental Wings	CoDiagnostix	☑	☑	In progress	☑	In progress	Codiagnostix Guide	☑	☑	☑	☑
Dentsply Sirona	Simplant	☑	☑	☑	☑	☑	Simplant SAFE guide (use NobelGuide*)	☑	–	☑	N/A
Digital Dental Service	DPS-pro	☑	☑	–	☑	–	DPS-pro	☑	☑	☑	☑
Dürr Dental	VistaSoft	☑	☑	–	☑	–	N/A				
Exocad	Exoplan	☑	☑	☑	☑	☑	Exoplan	In progress			
iCAT	LANDmarker	☑	☑	–	☑	–	Landmark Guide	☑	☑	☑	☑
Med3D	Implant 3D	☑	☑	–	–	–	N/A				
Media Lab	Implant 3D	☑	☑	In progress	☑	In progress	GuideDesign	☑	–	☑	☑
Morita	i-Dixel 2.0	☑	☑	In progress	☑	In progress	N/A				
MyRay	iRYS	☑	☑	–	–	–	N/A				
Nobel Biocare	DTX	☑	☑	☑	☑	☑	N/A				
NSequence	3Diagnosys	☑	☑	☑	☑	☑	RealGUIDE	☑	☑	☑	☑
	Maven Pro	☑	☑	–	–	–	NSequence Surgical Guide	☑	☑	☑	–
Planmeca	Planmeca Romexis	☑	☑	☑	☑	☑	3DDX Surgical Guide	☑	☑	☑	☑
ProDigiDent	ImplaStation	☑	☑	☑	☑	☑		–	–	–	–
SICAT	Galileos Implant	☑	☑	☑	☑	☑	SICAT surgical guides	☑	☑	☑	☑
	SICAT Implant	☑	☑	☑	☑	☑					
Swissmeda	SMOP	☑	☑	In progress	☑	In progress	SMOP	☑	☑	☑	☑
Zirkonzahn	Implant Planner	☑	☑	☑	☑	☑		–	–	–	–

* Use Replace RP 4.3 as reference implant and consider dimensional differences.

** Limitations apply depending on implant type, see instructions for use for details.

12112_3 2022-10

Neoss has provided information to the above specified software libraries if not stated otherwise. Neoss however is not the author, publisher, provider or controller of the software, the software provider or legal manufacturer of the surgical guides and as such does not support, guarantee the accuracy of, function of or compatibility of designs using this software or guides made by guide manufacturers for use with Neoss Guide Kit or any other Surgical Guide Kit. Neoss therefore has no liability to the Clinician or Patient for the use of this software or Surgical Guides.