



Your solution for digital
impressions *at the right price!*

NeoScan™ 1000



Intelligent Simplicity

Fast, accurate and *easy to use*

NeoScan™ 1000 is a fast, accurate and easy-to-use intraoral scanner that provides the possibility for a flexible workflow that saves you time and money.

- ✓ Easy to use
- ✓ Fast and accurate scanning
- ✓ Open system – implants, prosthetics and orthodontics
- ✓ Attractive price and no license fees

The simplified USB cable connection gives you the option to share the scanner between dental chairs in your practice. No separate power source is needed.

This easy-to-use scanner is what your practice needs. Best of all, the NeoScan 1000 comes at the right price and no license fees!



Digital dentistry *made easy*

When you purchase the NeoScan 1000, it includes the NeoPro software and NeoConnect cloud service with no additional license or cloud storage fees.

NeoPro™ software

NeoPro makes your digital impressions easy and fast. It's an open system with no software license fees. Amongst other it includes scan history management, smart reminders and free of charge software updates.



- ✓ Optimized scanning of metal/shining surfaces
- ✓ Automatic detection of margin lines
- ✓ Quadrant snapshots for aligner system
- ✓ Occlusal analysis and undercut check

NeoConnect® cloud service

NeoConnect is designed to improve workflow efficiency. It offers a seamless, efficient and secure process for data transfer and communication with your laboratory.



- ✓ No storage/cloud fees (up to 15 GB clinic, 5 GB lab)
- ✓ Add all the patient data you need
- ✓ Seamless process with your lab partners
- ✓ Full compliance with data privacy



The only digital impression machine you need!

You can use your NeoScan 1000 for your digital impression needs. It generates common STL and PLY file formats. These formats are industry standards and widely accepted in open systems.

Implant prosthetics

Open system for digital
implant impressions

ABUTMENTS,
BRIDGES,
AND BARS



Clear aligners

For design of 3D models
(STL/PLY format)

DIGITAL
ORTHODONTICS



SLEEP
APPLIANCES



CROWNS AND
BRIDGES WITH
PRE-PREPARATION
ON TOOTH LEVEL



Open system
for *C&B*
For all indications

Oral *sleep*
devices

For design of 3D models
(STL/PLY format)

Integrated
cloud service

For patient and lab
partner management



Better or equal to industry leading scanners




Michael Braian, PHD
Malmö, Sweden

Dr. Michael Braian's bullet chart shows the comparison of inter cylindrical precision on a dentated model. It shows that the NeoScan 1000 is better or equal to industry leading scanners on the market today. Watch his presentation from Neoss Integrate 2022 to learn more.



Webinar



Clinical case

Fully digital workflow using the Esthetic Healing Abutment with ScanPeg, intraoral scanning with the NeoScan 1000 and CAD design of the customized prosthetic restoration. Watch the procedure to learn more.



Clinical case



Marcus Dagnelid, DDS
Specialist in prosthodontics
Gothenburg, Sweden



Figure 1. A Neoss ProActive Edge implant in position 14 with an Esthetic Healing Abutment.



Figure 2. The ScanPeg is mounted on the Esthetic Healing Abutment. Intra oral scanning is performed with the NeoScan 1000.





Figure 3. The scan file is sent to the laboratory for fabrication of the final restoration.



Figure 4. Final restoration, customized abutment and crown, at time of delivery.

Digital Simplicity for you and your patients



Assortment & Warranties

Art. No.	Description
91000	NeoScan 1000 Intraoral Scanner <i>Includes Handpiece, 4 tips (2 x S, 2 x L), Holder Set, Manuals, Suitcase Extended warranties are available for purchase</i>
91012	NeoScan 1000 Intraoral Scanner <i>Includes 2 pcs Tip L TP102N</i>
91013	NeoScan 1000 Intraoral Scanner <i>Includes 2 pcs Tip S TP103N</i>
91022	Trolley for Intraoral Scanner

Digital libraries can be found under
<https://www.neoss.com/en-gb/cad-libraries>

Art. No.	Description
91024	NeoScan 1000 – Warranty extension – 1 year
91026	NeoScan 1000 – Warranty extension – 2 years
91027	NeoScan 1000 – Warranty extension – 3 years

UNITED KINGDOM & IRELAND OFFICE

Neoss UK & Ireland Ltd

Windsor House, Cornwall Road, Harrogate, HG1 2PW

T: +44 1423 817 733

E: info@neoss.com

E: NeoScan: neossupport.uk@neoss.com

For support, please visit our webpage:

neoss.com/neoscan1000

NeoScan customer support: +44 1423 533 991

Contact your local representative for product availability details. All products are not available in every market.

14825_2 EN 23-02 © Neoss Group, 2023. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss Group. Neoss and the N logo are registered trademarks of Neoss Group. NeoScan™, NeoPro™, and NeoConnect™ are trademarks of Neoss Group.