



Your in-lab  
ceramic solution

neoss® | NeoBase™ Abutments

Intelligent Simplicity

A complete range of titanium base abutments for local milling of ceramic restorations with angulated and straight screw channels.



### Full flexibility for individual needs

- Variety of post heights for additional retention area
- Variety of gingival margin heights for ideal implant placement
- All NeoBase™ abutments cater for both straight and angulated screw channels

### Angulated screw channels – improved handling and esthetics

- Favorable screw channel exit positions
- Screw-retained solution in the anterior for improved esthetics and occlusal function
- Easier posterior access



### Design freedom and production adaptability

- For use at angles up to 25° with 360° rotational freedom
  - Compact inner and outer geometry for maximum design freedom and compatibility with conventional milling machines
  - NeoBase™ post is designed for reliable retention and maximum ceramic thickness
-

Neoss provides a complete range of products to support your digital workflow and enable scanning of implant positions, design of prosthetics, and alternatives for local and central production.

3shape  exocad

### Design software and Libraries

Neoss provides design libraries for 3Shape and Exocad software to enable you to create your digital prosthetic design on NeoBase™ abutments.

Libraries together with the installation guidelines are available for download on [www.neoss-cadcam.de](http://www.neoss-cadcam.de).

### Design

A range of design options are available, including both single and multiple ceramic restorations, as well as the Esthetiline solution with matching gingival design and anatomical fit from healing abutments to final restorations.

### Milling

NeoBase™ mating geometry, including retentive features, is designed to maximize the options of compatible dental ceramic milling machines.

### Final result

NeoBase™ is the foundation for esthetically driven retrievable restorations, avoiding the need for intra-oral cementation. NeoBase™ is your local solution for ceramic restorations without compromise on design.



Neoss has been providing stress-free ceramic restorations for over 15 years with the original NeoLink® concept. The new range of NeoBase™ abutments caters for both straight and angulated screw channels, and is yet another example of our unique **Intelligent Simplicity** in action. Innovative and cost-effective solutions putting you, your team and your patients first.

# Assortment

Art. No.	Description
31367	NeoBase™ Mono G0.3 mm – H3.6 mm
31368	NeoBase™ Mono G1.5 mm – H3.6 mm
31369	NeoBase™ Mono G0.3 mm – H5.6 mm
31370	NeoBase™ Mono G1.5 mm – H5.6 mm
31371	NeoBase™ Multi G0.3 mm – H3.6 mm
<hr/>	
31372	NeoBase™ Mono G0.3 mm – H3.6 mm Ø3.25
31373	NeoBase™ Mono G0.3 mm – H5.6 mm Ø3.25
31374	NeoBase™ Multi G0.3 mm – H3.6 mm Ø3.25

Contact your local representative for information about availability in your market.

## Contact details

### Australia

Neoss Australia Pty. Ltd  
PO Box 404  
New Farm  
QLD 4005  
T +61 7 3216 0165  
F +61 7 3216 0135  
E info.au@neoss.com

### Austria

Neoss GmbH  
T +49 221 55405-322  
F +49 221 55405-522  
E info@neoss.de

### China

Shanghai Legend Life Technology Co., Ltd  
4th floor, Block 6  
1525 Minqiang Road  
Songjiang District  
Shanghai 201612  
P. R. China  
T +86-21-6765 7130-617  
F +86-21-6765 7131  
E info@legendlifetech.com

### Croatia

Medical Intertrade d.o.o.  
10431 Sveta Nedelja  
HRVATSKA  
T +385 1 33 74 020  
F +385 1 33 36 035  
E dental@medical-intertrade.hr

### Denmark

Neoss AB  
T +46 (0)31 88 12 80  
E info@neoss.se

### Germany

Neoss GmbH  
Im MediaPark 8  
DE-50670 Köln  
T +49 221 55405-322  
F +49 221 55405-522  
E info@neoss.de

### Italy

Neoss Italia S.r.l.  
Viale Certosa 138  
IT-20156 Milano  
T +39 02 92952 1  
F +39 02 92952 250  
E italia@neoss.info

### Japan

Neoss Japan Co.  
Y101 Building 5F  
No. 1-4 Kanda-Sudacho Chiyodaku  
Tokyo 101-0041  
T +81 3 5289 3511  
F +81 3 5289 3512  
E info@neoss.co.jp

### Middle East & Africa

Neoss Limited  
Unit 3737 Level 1  
Jewelry & Gemplex 3  
P.O. Box 5003315  
Dubai  
United Arab Emirates  
E info.mea@neoss.com

### Netherlands

Rotor, Medical Supplies  
Economiestraat 16  
NL-6135 KV Sittard  
T +31-46-481 12 03  
F +31-46-481 12 04  
E mail@rotor-medical.com

### New Zealand

Neoss Australia Pty. Ltd  
T +61 7 3216 0165  
F +61 7 3216 0135  
E info.au@neoss.com

### Norway

Neoss AB  
T +46 (0)31 88 12 80  
E info@neoss.se

### Republic of Ireland

Neoss Ltd  
T +44 1423 817-733  
F +44 1423 817-744  
E info@neoss.com

### Spain

Formación, Gestión Y Rendimiento, S.L  
C. Londres 54 E2B  
ES-08036 Barcelona  
T +34 931 408 770  
E info@fgrdental.com

### Sweden

Neoss AB  
Arvid Wallgrens backe 20  
SE-413 46 Göteborg  
T +46 (0)31 88 12 80  
E info@neoss.se

### Arc Solutions AB

Box 13011  
SE-250 13 Helsingborg  
T +46 (0)42 3017440  
E info@arcsolutions.se

### Switzerland

Flexident AG  
Amlehnstr. 22  
CH-6010 Kriens  
T +41413104020  
F +41413104025  
E info@flexident.ch

### Thailand

S.D.Tantavade Co., Ltd.  
47 Samakkee 20, Samakkee Road  
Thasai Muang  
Nonthanburi 11000  
T +66 2 952 4815 18  
F +66 2 952 4819  
E info@sdt1988.com

### Turkey

Dentmark Sağlık Ürünleri Anonim Şirketi  
Girne Mah. Girne Cad. No: 117  
Omega Deniz Plaza Kat: 2/12  
Maltepe  
Istanbul  
T +90 216 417 4040  
F +90 216 417 4041  
E info@neosshsi.com

### United Kingdom

Neoss Ltd  
Windsor House  
Cornwall Road  
Harrogate  
HG1 2PW  
T +44 1423 817-733  
F +44 1423 817-744  
E info@neoss.com

### United States

Neoss Inc.  
21860 Burbank Blvd.,  
Ste. 190 (North Entrance)  
Woodland Hills, CA 91367  
T +1 866 626-3677  
F +1 818 432-2640  
E contact.usa@neoss.com