



Fast, precise and *easy to use*

NeoScan™ 1000



Intelligent Simplicity

Fast scanning for accelerated results

The NeoScan™ 1000 is a fast, precise and easy-to-use intraoral scanner that provides the possibility for a flexible workflow. The simplified USB cable connection gives you the option to share the scanner between dental chairs in your practice. No separate power source is needed.

This easy-to-use scanner is what your practice needs. It supports full-touch screen and there is no need for any powder or pre-heating. The best of all is that the NeoScan 1000 comes at the right price!

- ✓ Easy to use
- ✓ Fast scanning
- ✓ High precision
- ✓ Attractive price



“The NeoScan 1000 is a *superfast, lightweight, and easy-to-use* scanner. I had the pleasure of being part of early testing and have used the scanner for several digital impression indications at my clinic with excellent results. Digital dentistry is in need of more *cost-efficient* solutions so that clinicians can use it to its full potential. The NeoScan 1000 has the potential to do just this.”

– Dr. Marcus Dagnelid
Specialist in prosthodontics | Gothenburg, Sweden



Your solution for *digital impressions*

Compact

The handle is ergonomically designed and is very easy to hold. Even with a compact size it achieves a 16 mm x 14 mm high-definition field of view.

Lightweight

With a handle weight of only 198 grams, like a smart phone, it will not cause discomfort even when held for a long time.

Vivid color

Thanks to high-resolution and true-color images, it allows you to distinguish tooth structure and soft tissues, so that margin lines and undercut areas can be easily identified.

Simplified connection

It is easy to share the scanner between dental chairs. NeoScan 1000 is powered directly by the computer, no need for a separate power supply.

Touch control

With a motion-capture feature and buttons on the handpiece, you can complete the entire process of scanning, optimizing, and exporting without touching the computer.

Attractive price

The best part is that this cutting-edge technology with easy-to-use features comes at a price level that is unbeatable.

Digital dentistry *made easy*

The NeoPro™ software enables you to produce restorations in your own practice or connect to industrial production. It makes your digital impression easy and fast. It's an open system and includes features such as scan history management, smart reminders, margin line tools, inspection and measurement.

Digital impression data sharing and order communication is designed to improve workflow efficiency. With the integrated cloud storage, NeoConnect™ offers a communication tool for data transfer and communication with your laboratory providing a seamless and efficient process.



Hear it from the *experts*



“I had the pleasure of being one of the beta testers for the NeoScan 1000 and I was positively surprised. The scanner is *superfast and lightweight, the data can be processed without any problems* and the first clinical results were convincing.”

– Dr. Fabia Flensburg
Dentist | Essen, Germany



“We have a large experience on intraoral scanners as we have had the opportunity to test the most important brands available on the market. The NeoScan 1000 is an *extremely fast and easy-to-use scanner*. On teeth and on implants it revealed a *great accuracy* in detecting all the details even in difficult clinical conditions.

The generated file is totally open so it is simple to share it with the clinical and technical components of the team. The NeoScan is an *essential device* for a full digital workflow in any modern dental office.”

– Dr. Giuseppe Luongo
Specialist in maxillofacial surgery | Rome, Italy



“As pioneers in the world of digital implantology, we are well-frequented beta testers of new intraoral scanners on the market. The NeoScan 1000 shows *very promising results in terms of resolution, speed, and feasibility*, and already has a permanent place in our digital workflow.”


– Dr. Nico Laube
Oral maxillofacial and plastic surgeon | Frankfurt, Germany



“I had the pleasure being part of the feasibility testing* for the NeoScan 1000 intraoral scanner since May 2021, the preliminary test results for *accuracy looks promising*.”

– Dr. Michael Braian
Dentist and dental technician | Malmö, Sweden

*Dr. Michael Braian's bullet chart indicating the comparison of inter cylindrical precision on a dentated model.



“The use of digital technology is an integral part of our clinic workflow. Central to this is the intraoral scanning process for which we require *precision, speed, ease of use and streamlined integration* with our laboratory providers. After extensive testing, I am confident that the NeoScan 1000 fulfills all these requirements and I have no hesitation in incorporating it into our clinic.”

– Dr. Kavitha N. Shah
Specialist in prosthodontics | London, United Kingdom



“True to the ethos of Neoss, the NeoScan scanner is an *extremely simple, compact and easy-to-use* piece of kit. It's great for dentists like myself wanting to introduce digital dentistry into their clinics. The set up was very straightforward and *the workflow between the lab and myself has been seamless*. Our patients have been impressed at how quick and efficient the scanner is at recording high-quality information about their mouths. It has improved patient engagement enormously and simplified treatment appointments.”

– Dr. Emilie Abraham
Specialist in prosthodontics | London, United Kingdom

Our offices

ARC SOLUTIONS

Arc Solutions
Box 13011
SE-250 13 Helsingborg
Sweden
T: +46 42 301 74 40
E: info@arcsolutions.se

AUSTRALIA & NEW ZEALAND

Neoss Australia Pty. Ltd
PO Box 404
New Farm
QLD 4005
Phone: +61 7 3216 0165
E: info.au@neoss.com

EUROPEAN DISTRIBUTION

Neoss AG
The Circle 23
8058 Zurich-Airport
Switzerland
E: distributors@neoss.com

GERMANY & AUSTRIA

Neoss GmbH
Im Mediapark 5b
50670 Köln
T: +49 221 96980 10
F: +49 221 96980 199
E: info@neoss.de

CHINA

Neoss China
7F, #303, Songhu Rd. (KIC Build #11)
Yangpu District
Shanghai, 200433
PR China
T: +86 21 80143135
E: china@neoss.com

ITALY

Neoss Italia S.r.l.
Viale Certosa 138
IT-20156 Milano
T: +39 02 92952 1
E: italia@neoss.info

MIDDLE EAST & AFRICA

Neoss Limited (DMCC BRANCH)
Unit 1505 Fortune Tower
Cluster C – JLT
Dubai
P.O. Box 64093
United Arab Emirates
E: info.mea@neoss.com

SWEDEN, DENMARK & NORWAY

Neoss AB
Arvid Wallgrens backe 20
SE-413 46 Göteborg
T: +46 31 88 12 80
E: info@neoss.se

UNITED KINGDOM & REPUBLIC OF IRELAND

Neoss Ltd
Windsor House
Cornwall Road
Harrogate
HG1 2PW
T: +44 1423 817 733
F: +44 1423 817 744
E: info@neoss.com

UNITED STATES

Neoss Inc.
21860 Burbank Blvd.,
Ste. 190 (North Entrance)
Woodland Hills, CA 91367
T: +1 866 626 3677
F: +1 818 432 2640
E: contact.usa@neoss.com

Our distributors

BAHRAIN

Cigalah Gulf Medical
Building 670, Block 309, Road 919
Salmaniya, Bahrain
T: +973 17232170
E: support@cgmed.com
maha@cgmed.com

CROATIA

Medical Intertrade d.o.o
10431 Sveta Nedelja
HRVATSKA
T: +385 1 33 74 020
F: +385 1 33 36 035
E: dental@medical-intertrade.hr

EGYPT

Dental Future
Building No 5
Mohamed Hassanein Heikal
St. Off Abbas El Akkad Street
1st Floor Flat 3
Cairo, Egypt
T: +20 109 6696449
E: mahmod@dentalfuturekw.com

IRAQ

Eblal Scientific Bureau
Al-Ameerat street
Almansour
Baghdad
Iraq
T: +964 780 411 0110
E: info@eblal.com

JAPAN

Neoss Japan Co.
Y101 Building 5F
No. 1-4 Kanda-Sudacho Chiyodaku
Tokyo 101-0041
T: +81 3 5289 3511
F: +81 3 5289 3512
E: info@neoss.co.jp

JORDAN

Advanced Intelligence for Medical Supplies & Equipment Ltd.
Al-Swefeyeh Paris str
Alhusainy Complex
Amman 11196
Jordan
T: +962 7 97389535
F: +962 6 5826797
E: fhaddam@advanceime.com

KINGDOM OF SAUDI ARABIA

Tamer Group
Mohammad Ibn Abdul Wahab Street
Al Ammariya District
P.O. Box 180
Jeddah, 21411
Kingdom of Saudi Arabia
T: +966 12 643 5600
F: +966 12 642 4141
E: alaa.wanly@tamergr.com

KUWAIT

Dental Future Company
Salmiya Salem el Mubark street
Suhaib complex office 10
Kuwait City
Kuwait
T: +965 51155997
E: Info@dentalfuturekw.com

LEBANON

Dentaltech Sarl
Brasilia Street
Hikmeh Building
First Basement
Baabda, Beirut – Lebanon
T: +961 5 950 707
F: +961 5 950 808

NETHERLANDS

Rotor, Medical Supplies B.V.
Rijksweg Zuid 51
NL-6071HV Swalmen
T: +31 46 481 12 02
F: +31 46 481 12 04
E: mail@rotor-medical.com

OMAN

Al Zahrawi Medical Supplies LLC.
Office No: 16 , Al Dar building opposite
Al Zubair automotive, street 61, Ghala
area.
Muscat, Oman
P.O. Box: 1627
T: +968 24587308
F: +968 24587487
E: info.oman@zahrawigroup.com

QATAR

Al Zahrawi Medical LLC
Building No. 36
Zone No. 39
Al Sadd Area
Street No.808
Doha, Qatar, P.O. Box 63065
T: +974 4431 872
F: +974 44365216
E: info.qatar@zahrawigroup.com

SPAIN

Formación, Gestión Y Rendimiento, S.L.
C. Nil Fabra 34-36.
Etslo 2 B
08012. Barcelona
Spain
T: +34 931 408 770
E: info@fgdental.com

SWITZERLAND

Flexident AG
Amlehstr. 22
CH-6010 Kriens
T: +41 41 3104020
E: info@flexident.ch

THAILAND

S.D. Tantavade Co., Ltd.
47 Samakkee 20, Samakkee Road Thasai
Muang
Nonthaburi 11000
Thailand
T: +66 2 952 4815 18
F: +66 2 952 4819
E: info@sdt1988.com

TURKEY

Dentmark Sağlık Ürünleri Anonim Şirketi
Cevizli Mah. Cevizli D-100 Güney
Yanyol Cad.
No:2 Kuriş Kule Kat:3 D:52
Kartal-Istanbul PO: 34865
T: +90 216 417 4040
E: info@dentmark.com.tr

UNITED ARAB EMIRATES

(DUBAI & NORTHERN EMIRATES)
Al Zahrawi Medical Supplies
1st Floor, Ford Middle East &
Africa Building
Al Wasl Office Complex, Sheikh Zayed
Road, Downtown
Dubai, United Arab Emirates
P.O. Box: 5973
T: +971 4 262 2728
F: +971 4 262 5506
E: info@zahrawigroup.com

UNITED ARAB EMIRATES (ABU DHABI)

Elwan Technical Supplies
Airport Street Abu Dhabi UAE
T: +971 26393292
E: info@elwandental.com

Dedicated support is available at neoscan.support.us@neoss.com

Contact your local representative for product availability details.
All products are not available in every market.

14825_1 EN 22-08 © Neoss Group, 2022. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss. Neoss, the N logo and NeoScan™ are trademarks of Neoss Group.