

NeoGen® Allograft Bone



The products pictured above are for illustration purposes only and are not to scale.

NeoGen® Allograft Mineralized Particulate 250–1000µm

CORTICAL/CANCELLOUS BLEND		CORTICAL
SIZE	ART. NO.	ART. NO.
0.25cc	64063	64069
0.5cc	64064	64070
1.0cc	64065	64071
2.0cc	64066	64072
2.5cc	64067	64073
5.0cc	64068	64074

NeoGen® Allograft is a mineralized dental particulate from 100% bone.

Product features:

- Available as a cortical or cortical and cancellous blend designed for both space maintenance and infiltration.
- An average particle size of 780µm for a multitude of applications and ease of handling.
- Terminal sterilization by low dose gamma radiation (10–6SAL).
- Available in six sizes, packaged in easy to use screw top jars.
- Maxillofacial and dental applications include, but are not limited to: Periodontal defects, socket preservation, sinus augmentation, and repair of bony defects.