



Off-the-shelf revolutionary *2-1 solution*  
Esthetic Healing Abutments with ScanPeg

Intelligent Simplicity

# A smarter, *more efficient patient-friendly* process

The Neoss® Esthetic Healing Abutments with ScanPeg is a combination of healing abutment and scan abutment in a revolutionary 2-in-1 solution. This is a patient-friendly, simplified and effortless procedure for accurate single tooth intraoral scanning and reduces the risk of soft-tissue damage with multiple healing abutment removals.

- ✔ Patient friendly
- ✔ Less invasive
- ✔ Easy handling
- ✔ Excellent esthetics



## Providing treatment with *fewer steps*

Neoss always strives at developing products that provide a treatment that is less invasive and can increase patient comfort. The Neoss Esthetic Healing Abutments with ScanPeg is part of the NeossONE solution, a prosthetic platform where the same abutment fits all Neoss implants no matter the implant design or diameter. Providing a complete range of anatomically shaped PEEK Healing Abutments with a dedicated Scan Abutment, for simplified and accurate intraoral scanning.

### Easier handling provides, faster improved procedure

- Minimized component exchange
- Streamlined intraoral scanning
- Less chair time, fewer treatment steps



### More patient friendly due to undisturbed healing

- Biological seal is maintained, tissue level is preserved
- Less invasive for increased patient comfort

### Providing excellent esthetics

- Shapes soft tissue during healing
- Perfectly matched with Neoss Esthetiline final abutments and Neoss customized abutments

## A step ahead

### 01 Healing abutment placement


The Esthetic Healing Abutment is placed during first or second stage surgery and then left until the final abutment is placed.



### 02 ScanPeg placement and scanning

An accurate intraoral scan is carried out at healing abutment level utilizing the ScanPeg. The ScanPeg comes co-packed with the Esthetic Healing Abutment and is easy to both attach and remove.

The intraoral digital impression can be taken whenever you choose. ScanPeg shape is included in Neoss CAD libraries.




### 03 Placement of the final restoration

Flexibility in when you take the digital impression allows design work to continue during the healing phase. The healing abutment is not removed until the final restoration is placed, eliminating one complete surgical intervention.

Contours of Esthetic Healing Abutments are included in Neoss CAD libraries to design the emergence profile from patients contour to final restoration.


**exocad 3shape**



## Assortment

Article no.	Description
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor
31362	Esthetic Healing Abutment with ScanPeg – Canine
31363	Esthetic Healing Abutment with ScanPeg – Premolar
31364	Esthetic Healing Abutment with ScanPeg – Molar
31365	ScanPeg – 5 pcs

Part of the Neoss Esthetiline range



## Hear it from the *experts*



“The NeoScan 1000 is *very fast, easy to use, and delivers extremely precise data*. Together with the Esthetic Healing Abutment with ScanPeg, the digital workflow is raised to a new level: *Greater patient comfort, more time efficiency and undisturbed healing for better esthetics and function.*”

– Dr. Fabia Flensberg  
Dentist | Essen, Germany



“The Esthetic Healing Abutments with ScanPeg are an excellent addition to the Neoss range. *They simplify the restorative process* and ensure you don't need to disrupt the soft tissue. The ability of the lab to copy the esthetic emergence profile makes fitting restorations simple and *reduces excessive pressure on the soft tissues.*”

– Dr. Deepesh Patel  
Specialist in prosthodontics | Cambridge, United Kingdom



“With the ScanPeg, a surgeon or prosthodontist can make a digital impression at any point in time without disturbing the soft-tissue healing. *It is more than simple, and absolutely intelligent.*”

– Dr. Jakob Zwaan  
Dentist | Calusco D'Adda, Italy



“The Esthetic Healing Abutment has completely transformed the way in which we can control the emergence healing of a recently placed implant with *minimal effort and maximum results*. This is the only fully integrated system made specifically for digital dentistry, the ScanPeg truly allows for *one abutment, one time concept*. This is the best implementation of a restorative system in dentistry right now.”

– Dr. Robert Ritter  
Dentist | Florida, United States



“The NeoScan 1000 is a *superfast, lightweight, and easy-to-use* scanner. I had the pleasure of being part of early testing and have used the scanner for several digital impression indications at my clinic with excellent results. Digital dentistry is in need of more *cost-efficient* solutions so that clinicians can use it to its full potential. The NeoScan 1000 has the potential to do just this.”

– Dr. Marcus Dagnelid  
Specialist in prosthodontics | Gothenburg, Sweden

## Your solution for *digital impressions*

The NeoScan™ 1000 is a fast, precise and easy-to-use intraoral scanner that provides the possibility for a flexible workflow. The simplified USB cable connection gives you the option to share the scanner between dental chairs in your practice. No separate power source is needed.

This easy-to-use scanner fits very well with the Neoss Esthetic Healing Abutments and is what your practice needs. The best of all is that the NeoScan 1000 comes at the right price!

- ✓ Easy to use
- ✓ Fast scanning
- ✓ High precision
- ✓ Attractive price



## Our offices

### ARC SOLUTIONS

Arc Solutions  
Box 13011  
SE-250 13 Helsingborg  
Sweden  
T: +46 42 301 74 40  
E: info@arcsolutions.se

### AUSTRALIA & NEW ZEALAND

Neoss Australia Pty. Ltd  
PO Box 404  
New Farm  
QLD 4005  
Phone: +61 7 3216 0165  
E: info.au@neoss.com

### CHINA

Neoss China  
7F, #303, Songhu Rd. (KIC Build #11)  
Yangpu District  
Shanghai, 200433  
PR China  
T: +86 21 80143135  
E: china@neoss.com

### DISTRIBUTION

Neoss AG  
The Circle 23  
8058 Zurich-Airport  
Switzerland  
E: distributors@neoss.com

### GERMANY & AUSTRIA

Neoss GmbH  
Im Mediapark 5b  
50670 Köln  
T: +49 221 96980 10  
F: +49 221 96980 199  
E: info@neoss.de

### ITALY

Neoss Italia S.r.l.  
Viale Certosa 138  
IT-20156 Milano  
T: +39 02 92952 1  
E: italia@neoss.info

### MIDDLE EAST & AFRICA Neoss Limited (DMCC BRANCH)

Unit 1505 Fortune Tower  
Cluster C – JLT  
Dubai  
P.O. Box 64093  
United Arab Emirates  
E: info.mea@neoss.com

### SWEDEN, DENMARK & NORWAY

Neoss AB  
Arvid Wallgrens backe 20  
SE-413 46 Göteborg  
T: +46 31 88 12 80  
E: info@neoss.se

### UNITED KINGDOM & REPUBLIC OF IRELAND

Neoss Ltd  
Windsor House  
Cornwall Road  
Harrogate  
HG1 2PW  
T: +44 1423 817 733  
F: +44 1423 817 744  
E: info@neoss.com

### UNITED STATES

Neoss Inc.  
21860 Burbank Blvd.,  
Ste. 190 (North Entrance)  
Woodland Hills, CA 91367  
T: +1 866 626 3677  
F: +1 818 432 2640  
E: contact.usa@neoss.com

## Our distributors

### BAHRAIN

Cigalah Gulf Medical  
Building 670, Block 309, Road 919  
Salmaniya, Bahrain  
T: +973 17232170  
E: support@cgmed.com  
maha@cgmed.com

### CROATIA

Medical Intertrade d.o.o  
10431 Sveta Nedelja  
HRVATSKA  
T: +385 1 33 74 020  
F: +385 1 33 36 035  
E: dental@medical-intertrade.hr

### EGYPT

Dental Future  
Building No 5  
Mohamed Hassanein Heikal  
St. Off Abbas El Akkad Street  
1st Floor Flat 3  
Cairo, Egypt  
T: +20 109 6696449  
E: mahmod@dentalfuturekw.com

### IRAQ

Eblal Scientific Bureau  
Al-Ameerat street  
Almansour  
Baghdad  
Iraq  
T: +964 780 411 0110  
E: info@eblal.com

### JAPAN

Neoss Japan Co.  
Y101 Building 5F  
No. 1-4 Kanda-Sudacho Chiyodaku  
Tokyo 101-0041  
T: +81 3 5289 3511  
F: +81 3 5289 3512  
E: info@neoss.co.jp

### JORDAN

Advanced Intelligence for Medical  
Supplies & Equipment Ltd.  
Al-Swefeyeh Paris str  
Alhusainy Complex  
Amman 11196  
Jordan  
T: +962 7 97389535  
F: +962 6 5826797  
E: f.khaddam@advanceime.com

### KINGDOM OF SAUDI ARABIA

Tamer Group  
Mohammad Ibn Abdul Wahab Street  
Al Ammariya District  
P.O. Box 180  
Jeddah, 21411  
Kingdom of Saudi Arabia  
T: +966 12 643 5600  
F: +966 12 642 4141  
E: alaa.wanly@tamergruop.com

### KUWAIT

Dental Future Company  
Salmiya Salem el Mubark street  
Suhaib complex office 10  
Kuwait City  
Kuwait  
T: +965 51155997  
E: Info@dentalfuturekw.com

### LEBANON

Dentaltech Sarl  
Brasilia Street  
Hikmeh Building  
First Basement  
Baabda, Beirut – Lebanon  
T: +961 5 950 707  
F: +961 5 950 808

### NETHERLANDS

Rotor, Medical Supplies B.V.  
Rijksweg Zuid 51  
NL-6071HV Swalmen  
T: +31 46 481 12 02  
F: +31 46 481 12 04  
E: mail@rotor-medical.com

### OMAN

Al Zahrawi Medical Supplies LLC.  
Office No: 16 , Al Dar building opposite  
Al Zubair automotive, street 61, Ghala  
area.  
Muscat, Oman  
P.O. Box: 1627  
T: +968 24587308  
F: +968 24587487  
E: info.oman@zahrawigroup.com

### QATAR

Al Zahrawi Medical LLC  
Building No. 36  
Zone No. 39  
Al Sadd Area  
Street No.808  
Doha, Qatar, P.O. Box 63065  
T: +974 4431 872  
F: +974 44365216  
E: info.qatar@zahrawigroup.com

### SPAIN

Formación, Gestión Y Rendimiento, S.L.  
C. Nil Fabra 34-36.  
Etslo 2 B  
08012. Barcelona  
Spain  
T: +34 931 408 770  
E: info@fgdental.com

### SWITZERLAND

Flexident AG  
Amlehstr. 22  
CH-6010 Kriens  
T: +41 41 3104020  
E: info@flexident.ch

### THAILAND

S.D. Tantavade Co., Ltd.  
47 Samakkee 20, Samakkee Road Thasai  
Muang  
Nonthaburi 11000  
Thailand  
T: +66 2 952 4815 18  
F: +66 2 952 4819  
E: info@sdt1988.com

### TURKEY

Dentmark Sağlık Ürünleri  
Anonim Şirketi  
Cevizli Mah. Cevizli D-100 Güney  
Yanyol Cad.  
No:2 Kuriş Kule Kat:3 D:52  
Kartal-Istanbul PO: 34865  
T: +90 216 417 4040  
E: info@dentmark.com.tr

### UNITED ARAB EMIRATES (ABU DHABI)

Elwan Technical Supplies  
Airport Street Abu Dhabi UAE  
T: +971 26393292  
E: info@elwandental.com

### UNITED ARAB EMIRATES (DUBAI & NORTHERN EMIRATES)

Al Zahrawi Medical Supplies  
1st Floor, Ford Middle East &  
Africa Building  
Al Wasl Office Complex, Sheikh Zayed  
Road, Downtown  
Dubai, United Arab Emirates  
P.O. Box: 5973  
T: +971 4 262 2728  
F: +971 4 262 5506  
E: info@zahrawigroup.com

Contact your local representative for product availability details.  
All products are not available in every market.

11885\_1 EN 22-09 © Neoss Group, 2022. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss. Neoss, the N logo and NeoScan™ are trademarks of Neoss Group.