

In-lab digital prosthetic

NeoBase® Abutments

- ✓ **Design flexibility** – Allows dental laboratories more freedom in manufacturing customized restorations.
- ✓ **Digital workflow** – Integration facilitate seamless design and anatomical fit from healing abutments to final restorations.
- ✓ **Compact geometry** – Offers reliable retention and maximize ceramic thickness.



Assortment

NeoBase® Abutments SSC (Straight Screw Channel)

All NeoBase Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	Abutment height	Gingival height	Diameter
31393	NeoBase Mono SSC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31394	NeoBase Mono SSC G1.5mm – H3.6mm	5.1mm	1.5mm	
31395	NeoBase Mono SSC G0.3mm – H5.6mm	5.9mm	0.3mm	
31396	NeoBase Mono SSC G1.5mm – H5.6mm	7.1mm	1.5mm	
31397	NeoBase Multi SSC G0.3mm – H3.6mm	3.9mm	0.3mm	
31398	NeoBase Mono SSC G0.3mm – H3.6mm, NP	3.9mm	0.3mm	3.8mm
31399	NeoBase Mono SSC G0.3mm – H5.6mm, NP	5.9mm		
31400	NeoBase Multi SSC G0.3mm – H3.6mm, NP	3.9mm		

Include Neo Abutment Screw NP and Neo Laboratory Screw NP

NeoBase® Abutments ASC (Angulated Screw Channel)

All NeoBase Abutments ASC include iGO Abutment Screw and iGO Laboratory Screw



Art. No.	Description	Abutment height	Gingival height	Diameter
31367	NeoBase Mono ASC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31368	NeoBase Mono ASC G1.5mm – H3.6mm	5.1mm	1.5mm	
31369	NeoBase Mono ASC G0.3mm – H5.6mm	5.9mm	0.3mm	
31370	NeoBase Mono ASC G1.5mm – H5.6mm	7.1mm	1.5mm	
31371	NeoBase Multi ASC G0.3mm – H3.6mm	3.9mm	0.3mm	
31372	NeoBase Mono ASC G0.3mm – H3.6mm, NP	3.9mm	0.3mm	3.8mm
31373	NeoBase Mono ASC G0.3mm – H5.6mm, NP	5.9mm		
31374	NeoBase Multi ASC G0.3mm – H3.6mm, NP	3.9mm		

Include iGO Abutment Screw NP and iGO Laboratory Screw NP

Design libraries available for 3Shape® and exocad® software. To download, scan the QR code.



Design libraries





ARC
A Neoss® Group Company



Neoss® Custom Prosthetics

– Powered by ARC

ARC milling center specializing for over 17 years in manufacturing custom-made implant prosthetics. Offering outstanding precision and reliable quality from single to full-arch restorations.*

*Angulated Screw Channel is available for all custom-made products.

Customized abutments

Available in titanium, cobalt chromium and zirconia



Custom Abutments SSC



Zi solutions on Ti-Base SSC/ASC



Custom Abutments ASC

Bars and bridges

Milled in titanium, cobalt chromium, or in additive cobalt chromium



Partial Bridges SSC/ASC



Full-Arch Bridges SSC/ASC



Additive Bridges SSC/ASC



Vertex Horix Bars



Vertex Bars



Preci classic Bar

	Custom Abutment			Bridge		Bar	
	Ti	CoCr	ZrO	Ti	CoCr	Ti	CoCr
Neoss SP	☑	☑	☑	☑	☑	☑	☑
Neoss NP	☑	☑	☑	☑	☑	☑	☑
Neoss Access				☑	☑	☑	☑

To place your order, visit www.neoss.com

12405_2 USA 2023-02 © Neoss Group, 2023. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss. Neoss, the N logo and NeoBase® are trademarks of Neoss Group. exocad® is a registered trademark of exocad GmbH and 3Shape® is a registered trademark of 3Shape A/S.