

*More than* just a TiBase

## NeoBase<sup>®</sup> Abutments

NeoBase; Neoss' efficient and flexible Ti-Base concept, empowers you with an unparalleled range of restorative options catering for both straight (SSC) and angulated screw channels (ASC). From single-tooth to fully edentulous indications, NeoBase covers the entire spectrum, delivering efficient solutions for both straightforward and demanding esthetic cases. NeoBase ensures seamless integration into your treatment plan.

- ✓ **Unmatched flexibility** – Choose from a wide range of restorative options, from single-tooth to fully edentulous indications, and achieve optimal outcomes for your patients.
- ✓ **Streamlined efficiency** – With our EasyDigital workflow, experience increased productivity in your lab without compromising on quality.
- ✓ **Original Neoss Implant connection** – Benefit from the trusted Neoss implant connection, ensuring compatibility, stability, and peace of mind by knowing you are covered by our warranty.



## Assortment


### NeoBase® Abutments SSC (Straight Screw Channel)

also referred to as TiBase

All NeoBase® Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	Abutment height	Gingival height	Diameter
31393	NeoBase Mono SSC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31394	NeoBase Mono SSC G1.5mm – H3.6mm	5.1mm	1.5mm	
31395	NeoBase Mono SSC G0.3mm – H5.6mm	5.9mm	0.3mm	
31396	NeoBase Mono SSC G1.5mm – H5.6mm	7.1mm	1.5mm	
31397	NeoBase Multi SSC G0.3mm – H3.6mm	3.9mm	0.3mm	3.8mm
31398	NeoBase Mono SSC G0.3mm – H3.6mm, NP	3.9mm	0.3mm	
31399	NeoBase Mono SSC G0.3mm – H5.6mm, NP	5.9mm		
31400	NeoBase Multi SSC G0.3mm – H3.6mm, NP	3.9mm		

*Include Neo Abutment Screw NP and Neo Laboratory Screw NP*



### NeoBase® Abutments ASC (Angulated Screw Channel)

also referred to as TiBase

All NeoBase® Abutments ASC include iGO Abutment Screw and iGO Laboratory Screw



Art. No.	Description	Abutment height	Gingival height	Diameter
31367	NeoBase Mono ASC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31368	NeoBase Mono ASC G1.5mm – H3.6mm	5.1mm	1.5mm	
31369	NeoBase Mono ASC G0.3mm – H5.6mm	5.9mm	0.3mm	
31370	NeoBase Mono ASC G1.5mm – H5.6mm	7.1mm	1.5mm	
31371	NeoBase Multi ASC G0.3mm – H3.6mm	3.9mm	0.3mm	3.8mm
31372	NeoBase Mono ASC G0.3mm – H3.6mm, NP	3.9mm	0.3mm	
31373	NeoBase Mono ASC G0.3mm – H5.6mm, NP	5.9mm		
31374	NeoBase Multi ASC G0.3mm – H3.6mm, NP	3.9mm		

*Include iGO Abutment Screw NP and iGO Laboratory Screw NP*

CAD libraries available for 3Shape® and exocad® software. To download, scan the QR code.



exocad  
3shape

