

*Optimized flexibility* with a wide range of angles and heights


## Access Abutments

The Access Abutment offers an axial, angulated, and straight extension to the implant, making it easier to work on and restore at the abutment level rather than directly on the implant.



- ✓ **Esthetic low profile** – Balances strength and esthetics with a unique internal connection requiring a minimum of interocclusal height.
- ✓ **Clinical screw fixed within the abutment** – Our patented design secures the clinical screw in the abutment eliminating the risk of dropping and losing the screw.
- ✓ **Prosthetic flexibility** – Option between implant and abutment level restorative protocol in the same case including the utilization of an Access Abutment with internal connection to alter platform height or angle.






## Access Abutments

Clinical screw fixed within the Multi-Unit Abutment

Art. No.	Description	A	B
<i>All Access Abutments include Abutment and a Healing Abutment</i>			
31238	Access Abutment 1.5 mm	1.5 mm	4 mm
31215	Access Abutment 3.0 mm	3 mm	
31216	Access Abutment 4.0 mm	4 mm	
31218	Access Abutment 6.0 mm	6 mm	
31219	Access Abutment 10° 2.6 mm	2.6 mm	4 mm
31243	Access Abutment 10° 4.6 mm	4.6 mm	
31221	Access Abutment 20° 2.6 mm	2.6 mm	
31244	Access Abutment 20° 4.6 mm	4.6 mm	
31223	Access Abutment 30° 2.9 mm	2.9 mm	
31245	Access Abutment 30° 4.9 mm	4.9 mm	
31225	Access Impression Coping 8 mm and Replica <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>		
31366	Access Impression Coping 13 mm and Replica		
90183	Scan Body Neoss® Access (IO 1B-A)		



Art. No.	Description	A	B	C
31226	Access Ti NeoLink® Multi	4.3 mm	0.3 mm	4.5 mm
31227	Access Gold NeoLink® Multi <i>31226 and 31227 include NeoLink®, Prosthetic Screw, Provisional/Laboratory Screw and a Plastic Coping</i>	4.7 mm	0.7 mm	4.5 mm
31334	Access Ball Abutment 1 mm			
90207	Access Equator Abutment Kit 0.5 mm			
90208	Access Equator Abutment Kit 1.5 mm			
31301	Access Provisional Ti Abutment Multi <i>Includes Provisional/Laboratory Screw and Straight Coping</i>	15.7 mm	0.7 mm	4.5 mm
31266	Access Burnout Abutment Multi – 2 pcs <i>Includes Straight Coping</i>			

