

More than just a TiBase

NeoBase[®] Abutments

NeoBase; Neoss' efficient and flexible Ti-Base concept, empowers you with an unparalleled range of restorative options catering for both straight (SSC) and angulated screw channels (ASC). From single-tooth to fully edentulous indications, NeoBase covers the entire spectrum, delivering efficient solutions for both straightforward and demanding esthetic cases. NeoBase ensures seamless integration into your treatment plan.

- ✓ **Unmatched flexibility** – Choose from a wide range of restorative options, from single-tooth to fully edentulous indications, and achieve optimal outcomes for your patients.
- ✓ **Streamlined efficiency** – With our EasyDigital workflow, experience increased productivity in your lab without compromising on quality.
- ✓ **Original Neoss Implant connection** – Benefit from the trusted Neoss implant connection, ensuring compatibility, stability, and peace of mind by knowing you are covered by our warranty.



Assortment


NeoBase® Abutments SSC (Straight Screw Channel)

also referred to as TiBase

All NeoBase® Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw

| Art. No. | Description | Abutment height | Gingival height | Diameter |
|----------|---------------------------------------|-----------------|-----------------|----------|
| 31393 | NeoBase Mono SSC G0.3mm – H3.6mm | 3.9mm | 0.3mm | 4.5mm |
| 31394 | NeoBase Mono SSC G1.5mm – H3.6mm | 5.1mm | 1.5mm | |
| 31395 | NeoBase Mono SSC G0.3mm – H5.6mm | 5.9mm | 0.3mm | |
| 31396 | NeoBase Mono SSC G1.5mm – H5.6mm | 7.1mm | 1.5mm | |
| 31397 | NeoBase Multi SSC G0.3mm – H3.6mm | 3.9mm | 0.3mm | 3.8mm |
| 31398 | NeoBase Mono SSC G0.3mm – H3.6mm, NP | 3.9mm | 0.3mm | |
| 31399 | NeoBase Mono SSC G0.3mm – H5.6mm, NP | 5.9mm | | |
| 31400 | NeoBase Multi SSC G0.3mm – H3.6mm, NP | 3.9mm | | |

Include Neo Abutment Screw NP and Neo Laboratory Screw NP



NeoBase® Abutments ASC (Angulated Screw Channel)

also referred to as TiBase

All NeoBase® Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw



| Art. No. | Description | Abutment height | Gingival height | Diameter |
|----------|---------------------------------------|-----------------|-----------------|----------|
| 31367 | NeoBase Mono ASC G0.3mm – H3.6mm | 3.9mm | 0.3mm | 4.5mm |
| 31368 | NeoBase Mono ASC G1.5mm – H3.6mm | 5.1mm | 1.5mm | |
| 31369 | NeoBase Mono ASC G0.3mm – H5.6mm | 5.9mm | 0.3mm | |
| 31370 | NeoBase Mono ASC G1.5mm – H5.6mm | 7.1mm | 1.5mm | |
| 31371 | NeoBase Multi ASC G0.3mm – H3.6mm | 3.9mm | 0.3mm | 3.8mm |
| 31372 | NeoBase Mono ASC G0.3mm – H3.6mm, NP | 3.9mm | 0.3mm | |
| 31373 | NeoBase Mono ASC G0.3mm – H5.6mm, NP | 5.9mm | | |
| 31374 | NeoBase Multi ASC G0.3mm – H3.6mm, NP | 3.9mm | | |

Include iGO Abutment Screw NP and iGO Laboratory Screw NP

CAD libraries available for 3Shape® and exocad® software. To download, scan the QR code.



CAD libraries

exocad
3shape

