


Experience *easy-to-use full arch flexibility* with predictable results using Neoss4+


Multi-Unit Abutments

The Multi-Unit Abutment is expertly crafted to restore function and esthetics in patients for full-arch restorations or partial bridge, especially when utilizing the well-established Neoss4+ treatment solution backed by clinical and scientific evidence. The abutment incorporates a clinical screw for ease of use, and our proprietary design enables tightening up to 32 Ncm, ensuring a secure fit.



- ✓ **Clinical screw fixed within the abutment** – Our patented design secures the clinical screw in the abutment eliminating the risk of dropping it.
- ✓ **NeossONE** – One prosthetic screw and one screwdriver for all Neoss products.
- ✓ **10° option and lower vertical height** – Only possible with Neoss thanks to our patented design.
- ✓ **Esthetic yellow anodization** – Achieve better esthetics for peri-implant soft tissue with yellow colored anodized titanium.





Angled Multi-Unit Abutments come delivered with a holder


Assortment

Multi-Unit Abutments

Clinical screw fixed within the Multi-Unit Abutment
Healing Cap and Replica needs to be ordered separately

Art. No.	Description	A	B	C
31411	Multi-Unit Abutment 1.0 mm	1mm	4.6 mm	
31412	Multi-Unit Abutment 2.0 mm	2mm		
31413	Multi-Unit Abutment 3.0 mm	3mm		
31414	Multi-Unit Abutment 4.0 mm	4mm		
31425	Multi-Unit Abutment 10° 2.0 mm	2mm	4.6 mm	
31426	Multi-Unit Abutment 10° 3.0 mm	3mm		
31427	Multi-Unit Abutment 10° 4.0 mm	4mm		
31429	Multi-Unit Abutment 17° 2.5 mm	2.5mm	4.6 mm	
31430	Multi-Unit Abutment 17° 3.5 mm	3.5mm		
31431	Multi-Unit Abutment 17° 4.5 mm	4.5mm		
31433	Multi-Unit Abutment 30° 3.0 mm	3mm	4.6 mm	
31434	Multi-Unit Abutment 30° 4.0 mm	4mm		
31435	Multi-Unit Abutment 30° 5.0 mm	5mm		
31440	Multi-Unit Healing Cap	5mm	5mm	
31442	Multi-Unit Impression Coping Open tray 6mm*			
31443	Multi-Unit Impression Coping Closed tray*			
31444	Multi-Unit Replica – 5 pcs			
31445	Multi-Unit Provisional Abutment Ti	15.7mm	0.7mm	5mm
31446	Multi-Unit NeoBase® SSC G 0.5mm – H 4.5mm	5mm	0.5mm	5.1mm
31447	Multi-Unit NeoBase® ASC G 0.5mm – H 4.5mm iGO	5mm	0.5mm	5.1mm
31450	Multi-Unit Scan Body			
91074	Analog for Printed Model Neoss® Multi-Unit			

*Not co-packed and should be ordered separately



Multi-Unit Abutments



Multi-Unit Healing Cap



Multi-Unit Impression Coping Open tray



Multi-Unit Impression Coping Closed tray



Multi-Unit Replica



Multi-Unit NeoBase® SSC



Multi-Unit Provisional Abutment Ti



Multi-Unit NeoBase® ASC



Multi-Unit Scan Body



Analog for Printed Model Neoss® Multi-Unit